

PrograMill PM7

PrograMill PM7 Softwareupdate V84.78 *Changelog*

New features and bug fixes:

- **new feature:** new titanium tool implemented in the tool library (available soon)
- **new feature:** implementation of the new IPS e.max ZirCAD Prime picture
- **new feature:** B- and X-axes are moving during magazine change to have more space
- **new feature:** more user friendly handling of Digital Denture process
- **new feature:** Y-axis stops if material changer flap is opened during disc changing process
- **new feature:** easy sorting of jobs in the job list with the help of arrow keys
- **new feature:** added color codes in the tool selection of the Tool Library
- **new feature:** disc is presented after milling (B-axis rotates 30 degrees to the left)
- **new feature:** function to clean the axes in the settings menu
- **fixed critical bug:** "Auto Mode" issue
- **fixed critical bug:** User Interface does not start
- **fixed bug:** missing animation in spindle cleaning process added
- **fixed bug:** improved animation for removing tool from the spindle