

**Safety Data Sheet**  
acc. to OSHA HCS

Printing date 06/23/2015

Version number 6

Reviewed on 06/23/2015


## 1 Identification

- **Product identifier**
- **Trade name:** **Special High Fusing White Ceramic Solder**
- **Relevant identified uses of the substance or mixture and uses advised against**  
No further relevant information available.
- **Application of the substance / the mixture**  
Manufacture of dental prosthesis  
Dental alloy
- **Details of the supplier of the safety data sheet**
- **Manufacturer/Supplier:**  
Ivoclar Vivadent Inc.  
175 Pineview Drive, Amherst, N.Y. 14228  
USA  
Tel. +1 800 533 6825  
Fax +1 716 691 2285  
  
Ivoclar Vivadent Inc.  
1-6600 Dixie Road  
Mississauga, Ontario  
L5T 2Y2  
Canada  
Phone: +1 905 670 8499  
Fax: +1 905 670 3102
- **Information department:** Quality Assurance / Regulatory Affairs
- **Emergency telephone number:**  
24 Hour Emergency Assistance:  
Emergency-Call USA - Infotrac: 1-800-535-5053  
Emergency-Call Canada - Canutec: 1-613-996-6666  
  
General SDS Assistance:  
US: 1-800-533-6825  
Canada: 1-800-263-8182

## 2 Hazard(s) identification

- **Classification of the substance or mixture**  
Ox. Sol. 2 H272 May intensify fire; oxidizer.

---

- **Label elements**
- **GHS label elements**  
The product is classified and labeled according to the Globally Harmonized System (GHS).
- **Hazard pictograms**  
  
GHS03
- **Signal word** Danger
- **Hazard statements**  
May intensify fire; oxidizer.
- **Precautionary statements**  
Keep away from heat.  
Take any precaution to avoid mixing with combustibles.

(Contd. on page 2)

## Safety Data Sheet acc. to OSHA HCS

Printing date 06/23/2015

Version number 6

Reviewed on 06/23/2015

**Trade name: Special High Fusing White Ceramic Solder**

(Contd. of page 1)

Wear protective gloves / eye protection / face protection.

Keep/Store away from clothing/combustible materials.

In case of fire: Use for extinction: CO<sub>2</sub>, powder or water spray.


Dispose of contents/container in accordance with local/regional/national/international regulations.

· **Additional information:**

Alloys do not require a label, providing they do not present a hazard to human health by inhalation, ingestion or contact with skin or to the aquatic environment in the form in which they are placed on the market.

· **Classification system:**

· **NFPA ratings (scale 0 - 4)**



Health = 1

Fire = 0

Reactivity = 0

· **HMIS-ratings (scale 0 - 4)**



Health = 1

Fire = 0

Reactivity = 0

· **Other hazards**

· **Results of PBT and vPvB assessment**

· **PBT:** Not applicable.

· **vPvB:** Not applicable.

### 3 Composition/information on ingredients

· **Chemical characterization: Mixtures**

· **Description:** Dental alloy

· **Dangerous components:**

7440-22-4	silver	25-50%
7440-05-3	palladium	10-<25%
7440-74-6	indium	1-<2.5%

### 4 First-aid measures

· **Description of first aid measures**

· **General information:** No special measures required.

· **After inhalation:**

Grinding dust: Supply fresh air; consult doctor in case of complaints.

In case of unconsciousness place patient stably in side position for transportation.

· **After skin contact:**

After contact with the molten product, cool rapidly with cold water.

If skin irritation continues, consult a doctor.

· **After eye contact:**

Grinding dust: Rinse opened eye for several minutes under running water. If symptoms persist, consult a doctor.

Mechanical effects only.

· **After swallowing:** Seek medical treatment.

· **Information for doctor:**

· **Most important symptoms and effects, both acute and delayed** No further relevant information available.

(Contd. on page 3)

## Safety Data Sheet acc. to OSHA HCS

Printing date 06/23/2015

Version number 6

Reviewed on 06/23/2015

**Trade name: Special High Fusing White Ceramic Solder**

(Contd. of page 2)

- **Indication of any immediate medical attention and special treatment needed**  
No further relevant information available.

### 5 Fire-fighting measures

- **Extinguishing media**
- **Suitable extinguishing agents:** Use fire fighting measures that suit the environment.
- **Special hazards arising from the substance or mixture** No further relevant information available.
- **Advice for firefighters**
- **Protective equipment:** Do not inhale explosion gases or combustion gases.

### 6 Accidental release measures

- **Personal precautions, protective equipment and emergency procedures** Not required.
- **Environmental precautions:**  
No special measures required.  
Do not allow to enter sewers/ surface or ground water.
- **Methods and material for containment and cleaning up:** Pick up mechanically.
- **Reference to other sections**  
See Section 7 for information on safe handling.  
See Section 8 for information on personal protection equipment.  
See Section 13 for disposal information.

### 7 Handling and storage

- **Handling:**
- **Precautions for safe handling**  
Extractors are required on all machines used for thermal processing or splinter removal processes.
- **Information about protection against explosions and fires:** No special measures required.
- **Conditions for safe storage, including any incompatibilities**
- **Storage:**
- **Requirements to be met by storerooms and receptacles:** No special requirements.
- **Information about storage in one common storage facility:** Store away from reducing agents.
- **Further information about storage conditions:** None.
- **Specific end use(s)** No further relevant information available.

### 8 Exposure controls/personal protection

- **Additional information about design of technical systems:** No further data; see item 7.
- **Control parameters**

· **Components with limit values that require monitoring at the workplace:**

**7440-22-4 silver**

PEL Long-term value: 0.01 mg/m<sup>3</sup>

REL Long-term value: 0.01 mg/m<sup>3</sup>

TLV Long-term value: 0.1 mg/m<sup>3</sup>  
metal: dust and fume

(Contd. on page 4)

US

## Safety Data Sheet

acc. to OSHA HCS

Printing date 06/23/2015

Version number 6

Reviewed on 06/23/2015

**Trade name: Special High Fusing White Ceramic Solder**

(Contd. of page 3)

### 7440-74-6 indium

REL	Long-term value: 0.1 mg/m <sup>3</sup> as In
TLV	Long-term value: 0.1 mg/m <sup>3</sup> as In

- **Additional information:** The lists that were valid during the creation were used as basis.
- **Exposure controls**
- **Personal protective equipment:**
- **General protective and hygienic measures:**  
Usual hygienic measures for dental practice.  
Wash hands before breaks and at the end of work.  
Keep away from foodstuffs, beverages and feed.
- **Breathing equipment:** Use respiratory protective device against the effects of fumes/dust/aerosol.
- **Protection of hands:**  
Protective gloves should always be worn during mechanical and thermal processing.  
After use of gloves apply skin-cleaning agents and skin cosmetics.
- **Material of gloves**  
Mechanical processing:  
Leather gloves  
Strong gloves  
Thermal processing:  
Heat protection gloves
- **Penetration time of glove material**  
The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.
- **Eye protection:**  
Always wear safety goggles during mechanical processing (grinding, sawing /cutting, drilling, milling).

## 9 Physical and chemical properties

- **Information on basic physical and chemical properties**
- **General Information**
- **Appearance:**

Form:	Solid
Color:	Grey
Odor:	Odorless
Odour threshold:	Not determined.
- **pH-value:** Not applicable.
- **Change in condition**

Melting point/Melting range:	1045 - 1105 °C (1913 - 2021 °F)
Boiling point/Boiling range:	Undetermined.
- **Flash point:** Not applicable.
- **Flammability (solid, gaseous):** Not determined.
- **Auto igniting:** Product is not selfigniting.
- **Danger of explosion:** Product does not present an explosion hazard.
- **Explosion limits:**

Lower:	Not determined.
Upper:	Not determined.

(Contd. on page 5)

## Safety Data Sheet

acc. to OSHA HCS

Printing date 06/23/2015

Version number 6

Reviewed on 06/23/2015

**Trade name: Special High Fusing White Ceramic Solder**

(Contd. of page 4)

· <b>Vapor pressure:</b>	Not applicable.
· <b>Density:</b>	Not determined.
· <b>Relative density</b>	Not determined.
· <b>Vapour density</b>	Not applicable.
· <b>Evaporation rate</b>	Not applicable.
· <b>Solubility in / Miscibility with Water:</b>	Insoluble.
· <b>Partition coefficient (n-octanol/water):</b>	Not determined.
· <b>Viscosity:</b>	
<b>Dynamic:</b>	Not applicable.
<b>Kinematic:</b>	Not applicable.
· <b>Other information</b>	No further relevant information available.

### 10 Stability and reactivity

- **Reactivity** No further relevant information available.
- **Chemical stability** Stable under normal handling and storage conditions.
- **Thermal decomposition / conditions to be avoided:** No decomposition if used according to specifications.
- **Possibility of hazardous reactions** Reacts with reducing agents.
- **Conditions to avoid** No further relevant information available.
- **Incompatible materials:** No further relevant information available.
- **Hazardous decomposition products:** None under normal conditions of storage and use.

### 11 Toxicological information

- **Information on toxicological effects**
- **Acute toxicity:**
- **Primary irritant effect:**
- **on the skin:** No irritant effect.
- **on the eye:** Mechanical effects only.
- **Sensitization:** No sensitizing effects known.
- **Additional toxicological information:**  
Fumes or dusts generated from cutting or grinding operations may cause respiratory irritation.
- **Carcinogenic categories**

- **NTP (National Toxicology Program)**

None of the ingredients is listed.

- **OSHA-Ca (Occupational Safety & Health Administration)**

None of the ingredients is listed.

### 12 Ecological information

- **Toxicity**
- **Aquatic toxicity:** No further relevant information available.
- **Persistence and degradability** No further relevant information available.
- **Behavior in environmental systems:**
- **Bioaccumulative potential** No further relevant information available.
- **Mobility in soil** No further relevant information available.

(Contd. on page 6)

## Safety Data Sheet

acc. to OSHA HCS

Printing date 06/23/2015

Version number 6

Reviewed on 06/23/2015

**Trade name: Special High Fusing White Ceramic Solder**

(Contd. of page 5)

- **Additional ecological information:**
- **General notes:**  
Water hazard class 1 (Self-assessment): slightly hazardous for water  
Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.
- **Results of PBT and vPvB assessment**
- **PBT:** Not applicable.
- **vPvB:** Not applicable.
- **Other adverse effects** No further relevant information available.

### 13 Disposal considerations

- **Waste treatment methods**
- **Recommendation:**  
Take to an approved landfill or a waste incineration plant, under conditions approved by the local authority.
- **Uncleaned packagings:**
- **Recommendation:** Disposal must be made according to official regulations.

### 14 Transport information

- |  |  |
|--|--|
| · <b>UN-Number</b>   |  |
| · <b>DOT</b>   | NA3077   |
| · <b>ADR, ADN, IMDG, IATA</b>  | Void   |
| · <b>UN proper shipping name</b>   |  |
| · <b>ADR, ADN, IMDG, IATA</b>  | Void   |
| · <b>Transport hazard class(es)</b>  |  |
| · <b>DOT</b>   |  |
| · <b>Class</b>   | 9 Miscellaneous dangerous substances and articles                              |
| · <b>Label</b>   | 9  |
| · <b>ADR, ADN, IMDG, IATA</b>  |  |
| · <b>Class</b>   | Void   |
| · <b>Packing group</b>   |  |
| · <b>DOT</b>   | III  |
| · <b>ADR, IMDG, IATA</b>   | Void   |
| · <b>Environmental hazards:</b>  |  |
| · <b>Marine pollutant:</b>   | No   |
| · <b>Special precautions for user</b>  | Not applicable.  |
| · <b>Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code</b> | Not applicable.  |
| · <b>Transport/Additional information:</b>                                       | Product is not classified as a dangerous good for transport (ADR, IMDG, IATA). |
| · <b>UN "Model Regulation":</b>  | -  |

US

(Contd. on page 7)

## Safety Data Sheet

acc. to OSHA HCS

Printing date 06/23/2015

Version number 6

Reviewed on 06/23/2015

**Trade name: Special High Fusing White Ceramic Solder**

(Contd. of page 6)

### 15 Regulatory information

· **Safety, health and environmental regulations/legislation specific for the substance or mixture**

· **Sara**

· **Section 355 (extremely hazardous substances):**

None of the ingredients is listed.

· **Section 313 (Specific toxic chemical listings):**

7440-22-4 silver

7440-50-8 copper

· **TSCA (Toxic Substances Control Act):**

7440-57-5 Gold

7440-22-4 silver

7440-50-8 copper

7440-06-4 platinum

7440-74-6 indium

7439-88-5 iridium

· **Proposition 65**

· **Chemicals known to cause cancer:**

None of the ingredients is listed.

· **Chemicals known to cause reproductive toxicity for females:**

None of the ingredients is listed.

· **Chemicals known to cause reproductive toxicity for males:**

None of the ingredients is listed.

· **Chemicals known to cause developmental toxicity:**

None of the ingredients is listed.

· **Carcinogenic categories**

· **EPA (Environmental Protection Agency)**

7440-22-4 silver

D

7440-50-8 copper

D

· **TLV (Threshold Limit Value established by ACGIH)**

None of the ingredients is listed.

· **NIOSH-Ca (National Institute for Occupational Safety and Health)**

None of the ingredients is listed.

· **GHS label elements**

The product is classified and labeled according to the Globally Harmonized System (GHS).

· **Hazard pictograms**



GHS03

· **Signal word** Danger

· **Hazard statements**

May intensify fire; oxidizer.

(Contd. on page 8)

US

## Safety Data Sheet acc. to OSHA HCS

Printing date 06/23/2015

Version number 6

Reviewed on 06/23/2015

**Trade name: *Special High Fusing White Ceramic Solder***

(Contd. of page 7)

- **Precautionary statements**

*Keep away from heat.*

*Take any precaution to avoid mixing with combustibles.*

*Wear protective gloves / eye protection / face protection.*

*Keep/Store away from clothing/combustible materials.*

*In case of fire: Use for extinction: CO<sub>2</sub>, powder or water spray.*

*Dispose of contents/container in accordance with local/regional/national/international regulations.*

- **National regulations:**

- **Other regulations, limitations and prohibitive regulations**

*The product is a medical device according to the Directive 93/42/EEC.*

*This product is classified as a medical device under US and Canadian regulations and has been reviewed by the US Food and Drug Administration and Health Canada.*

- **Chemical safety assessment:** *A Chemical Safety Assessment has not been carried out.*

### 16 Other information

*This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.*

- **Date of preparation / last revision** 06/23/2015 / 5

- **Abbreviations and acronyms:**

*ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)*

*IMDG: International Maritime Code for Dangerous Goods*

*DOT: US Department of Transportation*

*IATA: International Air Transport Association*

*ACGIH: American Conference of Governmental Industrial Hygienists*

*EINECS: European Inventory of Existing Commercial Chemical Substances*

*ELINCS: European List of Notified Chemical Substances*

*CAS: Chemical Abstracts Service (division of the American Chemical Society)*

*NFPA: National Fire Protection Association (USA)*

*HMIS: Hazardous Materials Identification System (USA)*

*Ox. Sol. 2: Oxidising Solids, Hazard Category 2*

US