

**Safety Data Sheet**  
acc. to OSHA HCS

Printing date 06/24/2015

Version number 4

Reviewed on 06/23/2015

**1 Identification**

- **Product identifier**
- **Trade name:** **Universal Solder PKF**
- **Relevant identified uses of the substance or mixture and uses advised against**  
No further relevant information available.
- **Application of the substance / the mixture**  
Manufacture of dental prosthesis  
Dental alloy
- **Details of the supplier of the safety data sheet**
- **Manufacturer/Supplier:**  
Ivoclar Vivadent Inc.  
175 Pineview Drive, Amherst, N.Y. 14228  
USA  
Tel. +1 800 533 6825  
Fax +1 716 691 2285  
  
Ivoclar Vivadent Inc.  
1-6600 Dixie Road  
Mississauga, Ontario  
L5T 2Y2  
Canada  
Phone: +1 905 670 8499  
Fax: +1 905 670 3102
- **Information department:** Quality Assurance / Regulatory Affairs
- **Emergency telephone number:**  
24 Hour Emergency Assistance:  
Emergency-Call USA - Infotrac: 1-800-535-5053  
Emergency-Call Canada - Canutec: 1-613-996-6666  
  
General SDS Assistance:  
US: 1-800-533-6825  
Canada: 1-800-263-8182

**2 Hazard(s) identification**

- **Classification of the substance or mixture**  
The product is not classified according to the Globally Harmonized System (GHS).

---

- **Label elements**
- **GHS label elements** Void
- **Hazard pictograms** Void
- **Signal word** Void
- **Hazard statements** Void
- **Additional information:**  
Alloys do not require a label, providing they do not present a hazard to human health by inhalation, ingestion or contact with skin or to the aquatic environment in the form in which they are placed on the market.
- **Classification system:**
- **NFPA ratings (scale 0 - 4)**



(Contd. on page 2)

# Safety Data Sheet

## acc. to OSHA HCS

Printing date 06/24/2015

Version number 4

Reviewed on 06/23/2015

Trade name: **Universal Solder PKF**

(Contd. of page 1)

· **HMIS-ratings (scale 0 - 4)**

HEALTH	0	Health = 0
FIRE	0	Fire = 0
REACTIVITY	0	Reactivity = 0

- **Other hazards**
- **Results of PBT and vPvB assessment**
- **PBT:** Not applicable.
- **vPvB:** Not applicable.

### 3 Composition/information on ingredients

- **Chemical characterization: Mixtures**
- **Description:** Dental alloy

· **Dangerous components:**

7440-22-4	silver	25-50%
7440-66-6	zinc powder -zinc dust (stabilized)	2.5-<10%
7440-06-4	platinum	2.5-<10%

### 4 First-aid measures

- **Description of first aid measures**
- **General information:** No special measures required.
- **After inhalation:**  
Grinding dust: Supply fresh air; consult doctor in case of complaints.  
In case of unconsciousness place patient stably in side position for transportation.
- **After skin contact:**  
After contact with the molten product, cool rapidly with cold water.  
If skin irritation continues, consult a doctor.
- **After eye contact:**  
Grinding dust: Rinse opened eye for several minutes under running water. If symptoms persist, consult a doctor.  
Mechanical effects only.
- **After swallowing:** Seek medical treatment.
- **Information for doctor:**
- **Most important symptoms and effects, both acute and delayed** No further relevant information available.
- **Indication of any immediate medical attention and special treatment needed**  
No further relevant information available.

### 5 Fire-fighting measures

- **Extinguishing media**
- **Suitable extinguishing agents:**  
CO<sub>2</sub>, extinguishing powder or water spray. Fight larger fires with water spray or alcohol resistant foam.
- **Special hazards arising from the substance or mixture** No further relevant information available.
- **Advice for firefighters**
- **Protective equipment:** Do not inhale explosion gases or combustion gases.

US

(Contd. on page 3)

## Safety Data Sheet acc. to OSHA HCS

Printing date 06/24/2015

Version number 4

Reviewed on 06/23/2015

Trade name: **Universal Solder PKF**

(Contd. of page 2)

### 6 Accidental release measures

- **Personal precautions, protective equipment and emergency procedures** Not required.
- **Environmental precautions:**  
No special measures required.  
Do not allow to enter sewers/ surface or ground water.
- **Methods and material for containment and cleaning up:** Pick up mechanically.
- **Reference to other sections**  
See Section 7 for information on safe handling.  
See Section 8 for information on personal protection equipment.  
See Section 13 for disposal information.

### 7 Handling and storage

- **Handling:**
- **Precautions for safe handling**  
Extractors are required on all machines used for thermal processing or splinter removal processes.
- **Information about protection against explosions and fires:** No special measures required.
- **Conditions for safe storage, including any incompatibilities**
- **Storage:**
- **Requirements to be met by storerooms and receptacles:** No special requirements.
- **Information about storage in one common storage facility:** Store away from flammable substances.
- **Further information about storage conditions:** None.
- **Specific end use(s)** No further relevant information available.

### 8 Exposure controls/personal protection

- **Additional information about design of technical systems:** No further data; see item 7.
- **Control parameters**

· **Components with limit values that require monitoring at the workplace:**

**7440-22-4 silver**

PEL Long-term value: 0.01 mg/m<sup>3</sup>

REL Long-term value: 0.01 mg/m<sup>3</sup>

TLV Long-term value: 0.1 mg/m<sup>3</sup>  
metal: dust and fume

**7440-06-4 platinum**

PEL Long-term value: 0.002 mg/m<sup>3</sup>  
soluble salts, as Pt

REL Long-term value: 1\* 0.002\*\* mg/m<sup>3</sup>  
\*Metal \*\*soluble salts, as Pt

TLV Long-term value: 1\* 0.002\*\* mg/m<sup>3</sup>  
\*metal \*\*soluble salts, as Pt

- **Additional information:** The lists that were valid during the creation were used as basis.
- **Exposure controls**
- **Personal protective equipment:**
- **General protective and hygienic measures:**  
Usual hygienic measures for dental practice.  
Wash hands before breaks and at the end of work.  
Keep away from foodstuffs, beverages and feed.

(Contd. on page 4)

## Safety Data Sheet

### acc. to OSHA HCS

Printing date 06/24/2015

Version number 4

Reviewed on 06/23/2015

**Trade name: Universal Solder PKF**

(Contd. of page 3)

- **Breathing equipment:** Use respiratory protective device against the effects of fumes/dust/aerosol.
- **Protection of hands:**  
Protective gloves should always be worn during mechanical and thermal processing.  
After use of gloves apply skin-cleaning agents and skin cosmetics.
- **Material of gloves**  
Mechanical processing:  
Leather gloves  
Strong gloves  
Thermal processing:  
Heat protection gloves
- **Penetration time of glove material**  
The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.
- **Eye protection:**  
Always wear safety goggles during mechanical processing (grinding, sawing /cutting, drilling, milling).

### 9 Physical and chemical properties

- |   |   |
|---|---|
| · <b>Information on basic physical and chemical properties</b>    |   |
| · <b>General Information</b>                                      |   |
| · <b>Appearance:</b>  |   |
| Form:   | Solid   |
| Color:  | Grey  |
| · <b>Odor:</b>  | Odorless                                      |
| · <b>Odour threshold:</b>   | Not determined.                               |
| · <b>pH-value:</b>  | Not applicable.                               |
| · <b>Change in condition</b>                                      |   |
| Melting point/Melting range:                                      | 1100 - 1165 °C (2012 - 2129 °F)               |
| Boiling point/Boiling range:                                      | Undetermined.                                 |
| · <b>Flash point:</b>   | Not applicable.                               |
| · <b>Flammability (solid, gaseous):</b>                           | Not determined.                               |
| · <b>Auto igniting:</b>   | Product is not selfigniting.                  |
| · <b>Danger of explosion:</b>                                     | Product does not present an explosion hazard. |
| · <b>Explosion limits:</b>  |   |
| Lower:  | Not determined.                               |
| Upper:  | Not determined.                               |
| · <b>Vapor pressure:</b>  | Not applicable.                               |
| · <b>Density:</b>   | Not determined.                               |
| · <b>Relative density</b>   | Not determined.                               |
| · <b>Vapour density</b>   | Not applicable.                               |
| · <b>Evaporation rate</b>   | Not applicable.                               |
| · <b>Solubility in / Miscibility with</b>                         |   |
| Water:  | Insoluble.                                    |
| · <b>Partition coefficient (n-octanol/water):</b> Not determined. |   |
| · <b>Viscosity:</b>   |   |
| Dynamic:  | Not applicable.                               |
| Kinematic:  | Not applicable.                               |

(Contd. on page 5)

US

## Safety Data Sheet acc. to OSHA HCS

Printing date 06/24/2015

Version number 4

Reviewed on 06/23/2015

**Trade name: Universal Solder PKF**

(Contd. of page 4)

· **Other information**

No further relevant information available.

### 10 Stability and reactivity

- **Reactivity** No further relevant information available.
- **Chemical stability** Stable under normal handling and storage conditions.
- **Thermal decomposition / conditions to be avoided:** No decomposition if used according to specifications.
- **Possibility of hazardous reactions** Reacts with reducing agents.
- **Conditions to avoid** No further relevant information available.
- **Incompatible materials:** No further relevant information available.
- **Hazardous decomposition products:** None under normal conditions of storage and use.

### 11 Toxicological information

- **Information on toxicological effects**
- **Acute toxicity:**
- **Primary irritant effect:**
- **on the skin:** No irritant effect.
- **on the eye:** Mechanical effects only.
- **Sensitization:** No sensitizing effects known.
- **Additional toxicological information:**  
Fumes or dusts generated from cutting or grinding operations may cause respiratory irritation.
- **Carcinogenic categories**

· **NTP (National Toxicology Program)**

None of the ingredients is listed.

· **OSHA-Ca (Occupational Safety & Health Administration)**

None of the ingredients is listed.

### 12 Ecological information

- **Toxicity**
- **Aquatic toxicity:** No further relevant information available.
- **Persistence and degradability** No further relevant information available.
- **Behavior in environmental systems:**
- **Bioaccumulative potential** No further relevant information available.
- **Mobility in soil** No further relevant information available.
- **Additional ecological information:**
- **General notes:**  
Water hazard class 2 (Self-assessment): hazardous for water  
Do not allow product to reach ground water, water course or sewage system.  
Danger to drinking water if even small quantities leak into the ground.
- **Results of PBT and vPvB assessment**
- **PBT:** Not applicable.
- **vPvB:** Not applicable.
- **Other adverse effects** No further relevant information available.

US

(Contd. on page 6)

## Safety Data Sheet

acc. to OSHA HCS

Printing date 06/24/2015

Version number 4

Reviewed on 06/23/2015

Trade name: **Universal Solder PKF**

(Contd. of page 5)

### 13 Disposal considerations

- **Waste treatment methods**
- **Recommendation:**  
Take to an approved landfill or a waste incineration plant, under conditions approved by the local authority.
- **Uncleaned packagings:**
- **Recommendation:** Disposal must be made according to official regulations.

### 14 Transport information

- |  |  |
|--|--|
| · <b>UN-Number</b>   |  |
| · <b>DOT</b>   | NA3077   |
| · <b>ADR, ADN, IMDG, IATA</b>  | Void   |
| · <b>UN proper shipping name</b>   |  |
| · <b>DOT</b>   | Oxidizing solid, n.o.s.  |
| · <b>ADR</b>   | 1479 Oxidizing solid, n.o.s., ENVIRONMENTALLY HAZARDOUS                        |
| · <b>ADN, IMDG, IATA</b>   | Void   |
| · <b>Transport hazard class(es)</b>  |  |
| · <b>DOT</b>   |  |
| · <b>Class</b>   | 9 Miscellaneous dangerous substances and articles                              |
| · <b>Label</b>   | 9  |
| · <b>ADR, ADN, IMDG, IATA</b>  |  |
| · <b>Class</b>   | Void   |
| · <b>Packing group</b>   |  |
| · <b>DOT</b>   | III  |
| · <b>ADR, IMDG, IATA</b>   | Void   |
| · <b>Environmental hazards:</b>  |  |
| · <b>Marine pollutant:</b>   | No   |
| · <b>Special precautions for user</b>  | Not applicable.  |
| · <b>Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code</b> | Not applicable.  |
| · <b>Transport/Additional information:</b>                                       | Product is not classified as a dangerous good for transport (ADR, IMDG, IATA). |
| · <b>UN "Model Regulation":</b>  | UN1479, Oxidizing solid, n.o.s., ENVIRONMENTALLY HAZARDOUS, 5.1, II            |

### 15 Regulatory information

- **Safety, health and environmental regulations/legislation specific for the substance or mixture**
- **Sara**

· **Section 355 (extremely hazardous substances):**

None of the ingredients is listed.

· **Section 313 (Specific toxic chemical listings):**

7440-22-4 | silver

(Contd. on page 7)

## Safety Data Sheet

acc. to OSHA HCS

Printing date 06/24/2015

Version number 4

Reviewed on 06/23/2015

**Trade name: Universal Solder PKF**

(Contd. of page 6)

7440-50-8 copper

· **TSCA (Toxic Substances Control Act):**

7440-57-5 Gold

7440-22-4 silver

7440-50-8 copper

7440-06-4 platinum

7440-74-6 indium

7439-88-5 iridium

· **Proposition 65**

· **Chemicals known to cause cancer:**

None of the ingredients is listed.

· **Chemicals known to cause reproductive toxicity for females:**

None of the ingredients is listed.

· **Chemicals known to cause reproductive toxicity for males:**

None of the ingredients is listed.

· **Chemicals known to cause developmental toxicity:**

None of the ingredients is listed.

· **Carcinogenic categories**

· **EPA (Environmental Protection Agency)**

7440-22-4 silver

D

7440-50-8 copper

D

· **TLV (Threshold Limit Value established by ACGIH)**

None of the ingredients is listed.

· **NIOSH-Ca (National Institute for Occupational Safety and Health)**

None of the ingredients is listed.

· **GHS label elements** Void

· **Hazard pictograms** Void

· **Signal word** Void

· **Hazard statements** Void

· **National regulations:**

· **Other regulations, limitations and prohibitive regulations**

The product is a medical device according to the Directive 93/42/EEC.

This product is classified as a medical device under US and Canadian regulations and has been reviewed by the US Food and Drug Administration and Health Canada.

· **Chemical safety assessment:** A Chemical Safety Assessment has not been carried out.

### 16 Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

· **Date of preparation / last revision** 06/24/2015 / 3

· **Abbreviations and acronyms:**

ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

DOT: US Department of Transportation

IATA: International Air Transport Association

ACGIH: American Conference of Governmental Industrial Hygienists

(Contd. on page 8)

**Safety Data Sheet**  
**acc. to OSHA HCS**

Printing date 06/24/2015

Version number 4

Reviewed on 06/23/2015

**Trade name: *Universal Solder PKF***

*EINECS: European Inventory of Existing Commercial Chemical Substances*  
*ELINCS: European List of Notified Chemical Substances*  
*CAS: Chemical Abstracts Service (division of the American Chemical Society)*  
*NFPA: National Fire Protection Association (USA)*  
*HMIS: Hazardous Materials Identification System (USA)*

(Contd. of page 7)

US