


**1 Identification**

- **Product identifier**
- **Trade name:** **Capricorn**
- **Relevant identified uses of the substance or mixture and uses advised against**  
No further relevant information available.
- **Application of the substance / the mixture** Manufacture of dental prothesis
- **Details of the supplier of the safety data sheet**
- **Manufacturer/Supplier:**  
Ivoclar Vivadent Inc.  
175 Pineview Drive, Amherst, N.Y. 14228  
USA  
Tel. +1 800 533 6825  
Fax +1 716 691 2285
  
- Ivoclar Vivadent Inc.  
1-6600 Dixie Road  
Mississauga, Ontario  
L5T 2Y2  
Canada  
Phone: +1 905 670 8499  
Fax: +1 905 670 3102
- **Information department:** Quality Assurance / Regulatory Affairs
- **Emergency telephone number:**  
24 Hour Emergency Assistance:  
Emergency-Call USA - Infotrac: 1-800-535-5053  
Emergency-Call Canada - Canutec: 1-613-996-6666
  
- General SDS Assistance:  
US: 1-800-533-6825  
Canada: 1-800-263-8182


**2 Hazard(s) identification**

- **Classification of the substance or mixture**  
Ox. Sol. 2 H272 May intensify fire; oxidizer.  
Skin Corr. 1C H314 Causes severe skin burns and eye damage.  
Eye Dam. 1 H318 Causes serious eye damage.

---

  - **Label elements**
  - **GHS label elements**  
The product is classified and labeled according to the Globally Harmonized System (GHS).
  - **Hazard pictograms**
- 

GHS03



GHS05
- **Signal word** Danger
  - **Hazard-determining components of labeling:**  
gallium
  - **Hazard statements**  
May intensify fire; oxidizer.

## Safety Data Sheet

acc. to OSHA HCS

Printing date 06/23/2015

Version number 5

Reviewed on 05/07/2015

**Trade name: Capricorn**

(Contd. of page 1)

*Causes severe skin burns and eye damage.*

**· Precautionary statements**

*Keep away from heat.*

*Take any precaution to avoid mixing with combustibles.*

*Wear protective gloves / eye protection / face protection.*

*Wear eye protection / face protection.*

*Do not breathe dusts or mists.*

*Keep/Store away from clothing/combustible materials.*

*Wash thoroughly after handling.*

*If on skin (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower.*

*IF INHALED: Remove person to fresh air and keep comfortable for breathing.*

*If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do.*

*Continue rinsing.*

*Immediately call a poison center/doctor.*

*Specific treatment (see on this label).*

*If swallowed: Rinse mouth. Do NOT induce vomiting.*

*Wash contaminated clothing before reuse.*

*In case of fire: Use for extinction: CO<sub>2</sub>, powder or water spray.*

*Store locked up.*

*Dispose of contents/container in accordance with local/regional/national/international regulations.*

**· Additional information:**

*Alloys do not require a label, providing they do not present a hazard to human health by inhalation, ingestion or contact with skin or to the aquatic environment in the form in which they are placed on the market.*

**· Classification system:**

**· NFPA ratings (scale 0 - 4)**



Health = 1

Fire = 0

Reactivity = 0

**· HMIS-ratings (scale 0 - 4)**



Health = 1

Fire = 0

Reactivity = 0

**· Other hazards**

**· Results of PBT and vPvB assessment**

**· PBT:** Not applicable.

**· vPvB:** Not applicable.

### 3 Composition/information on ingredients

**· Chemical characterization: Mixtures**

**· Description:** Dental alloy

**· Dangerous components:**

7440-05-3	palladium	50-100%
7440-74-6	indium	2.5-<10%
7440-55-3	gallium	2.5-<10%
7440-22-4	silver	2.5-<10%

(Contd. on page 3)

## Safety Data Sheet acc. to OSHA HCS

Printing date 06/23/2015

Version number 5

Reviewed on 05/07/2015

Trade name: **Capricorn**

(Contd. of page 2)

### 4 First-aid measures

- **Description of first aid measures**
- **General information:** No special measures required.
- **After inhalation:**  
Grinding dust: Supply fresh air; consult doctor in case of complaints.  
In case of unconsciousness place patient stably in side position for transportation.
- **After skin contact:** After contact with the molten product, cool rapidly with cold water.
- **After eye contact:**  
Grinding dust: Rinse opened eye for several minutes under running water. If symptoms persist, consult a doctor.  
Mechanical effects only.
- **After swallowing:** Seek medical treatment.
- **Information for doctor:**
- **Most important symptoms and effects, both acute and delayed** No further relevant information available.
- **Indication of any immediate medical attention and special treatment needed**  
No further relevant information available.

### 5 Fire-fighting measures

- **Extinguishing media**
- **Suitable extinguishing agents:** The product is not flammable.
- **Special hazards arising from the substance or mixture** No further relevant information available.
- **Advice for firefighters**
- **Protective equipment:** Do not inhale explosion gases or combustion gases.

### 6 Accidental release measures

- **Personal precautions, protective equipment and emergency procedures** Not required.
- **Environmental precautions:** No special measures required.
- **Methods and material for containment and cleaning up:** Pick up mechanically.
- **Reference to other sections**  
See Section 7 for information on safe handling.  
See Section 8 for information on personal protection equipment.  
See Section 13 for disposal information.

### 7 Handling and storage

- **Handling:**
- **Precautions for safe handling**  
Extractors are required on all machines used for thermal processing or splinter removal processes.
- **Information about protection against explosions and fires:** No special measures required.
- **Conditions for safe storage, including any incompatibilities**
- **Storage:**
- **Requirements to be met by storerooms and receptacles:** No special requirements.
- **Information about storage in one common storage facility:** Not required.
- **Further information about storage conditions:** None.
- **Specific end use(s)** No further relevant information available.

US

(Contd. on page 4)

# Safety Data Sheet

## acc. to OSHA HCS

Printing date 06/23/2015

Version number 5

Reviewed on 05/07/2015

Trade name: **Capricorn**

(Contd. of page 3)

### 8 Exposure controls/personal protection

- **Additional information about design of technical systems:** No further data; see item 7.
- **Control parameters**

· **Components with limit values that require monitoring at the workplace:**

#### 7440-74-6 indium

REL	Long-term value: 0.1 mg/m <sup>3</sup> as In
TLV	Long-term value: 0.1 mg/m <sup>3</sup> as In

#### 7440-22-4 silver

PEL	Long-term value: 0.01 mg/m <sup>3</sup>
REL	Long-term value: 0.01 mg/m <sup>3</sup>
TLV	Long-term value: 0.1 mg/m <sup>3</sup> metal: dust and fume

- **Additional information:** The lists that were valid during the creation were used as basis.
- **Exposure controls**
- **Personal protective equipment:**
- **General protective and hygienic measures:** Usual hygienic measures for dental practice.
- **Breathing equipment:** Use respiratory protective device against the effects of fumes/dust/aerosol.
- **Protection of hands:**  
Protective gloves should always be worn during mechanical and thermal processing.  
After use of gloves apply skin-cleaning agents and skin cosmetics.
- **Material of gloves**  
Mechanical processing:  
Leather gloves  
Strong gloves  
Thermal processing:  
Heat protection gloves
- **Penetration time of glove material**  
The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.
- **Eye protection:**  
Always wear safety goggles during mechanical processing (grinding, sawing /cutting, drilling, milling).

### 9 Physical and chemical properties

· **Information on basic physical and chemical properties**

· **General Information**

· **Appearance:**

<b>Form:</b>	Solid
<b>Color:</b>	White
<b>Odor:</b>	Odorless
<b>Odour threshold:</b>	Not determined.

· **pH-value:** Not applicable.

· **Change in condition**

<b>Melting point/Melting range:</b>	1170–1335 °C (2138–2435 °F)
<b>Boiling point/Boiling range:</b>	Undetermined.

(Contd. on page 5)

US

## Safety Data Sheet

acc. to OSHA HCS

Printing date 06/23/2015

Version number 5

Reviewed on 05/07/2015

Trade name: **Capricorn**

(Contd. of page 4)

· <b>Flash point:</b>	Not applicable.
· <b>Flammability (solid, gaseous):</b>	Not determined.
· <b>Auto igniting:</b>	Product is not selfigniting.
· <b>Danger of explosion:</b>	Product does not present an explosion hazard.
· <b>Explosion limits:</b>	
<b>Lower:</b>	Not determined.
<b>Upper:</b>	Not determined.
· <b>Vapor pressure:</b>	Not applicable.
· <b>Density at 20 °C (68 °F):</b>	11 g/cm <sup>3</sup> (91.795 lbs/gal)
· <b>Relative density</b>	Not determined.
· <b>Vapour density</b>	Not applicable.
· <b>Evaporation rate</b>	Not applicable.
· <b>Solubility in / Miscibility with Water:</b>	Insoluble.
· <b>Partition coefficient (n-octanol/water):</b>	Not determined.
· <b>Viscosity:</b>	
<b>Dynamic:</b>	Not applicable.
<b>Kinematic:</b>	Not applicable.
· <b>Solvent content:</b>	
<b>Solids content:</b>	100 %
· <b>Other information</b>	No further relevant information available.

### 10 Stability and reactivity

- **Reactivity** No further relevant information available.
- **Chemical stability** Stable under normal handling and storage conditions.
- **Thermal decomposition / conditions to be avoided:** No decomposition if used according to specifications.
- **Possibility of hazardous reactions** No dangerous reactions known.
- **Conditions to avoid** No further relevant information available.
- **Incompatible materials:** No further relevant information available.
- **Hazardous decomposition products:** None under normal conditions of storage and use.

### 11 Toxicological information

- **Information on toxicological effects**
- **Acute toxicity:**
- **Primary irritant effect:**
- **on the skin:** No irritant effect.
- **on the eye:** No irritating effect.
- **Sensitization:** No sensitizing effects known.
- **Additional toxicological information:**  
Fumes or dusts generated from cutting or grinding operations may cause respiratory irritation.
- **Carcinogenic categories**

#### · **NTP (National Toxicology Program)**

None of the ingredients is listed.

(Contd. on page 6)

US

## Safety Data Sheet

acc. to OSHA HCS

Printing date 06/23/2015

Version number 5

Reviewed on 05/07/2015

Trade name: **Capricorn**

(Contd. of page 5)

· **OSHA-Ca (Occupational Safety & Health Administration)**

None of the ingredients is listed.

### 12 Ecological information

- **Toxicity**
- **Aquatic toxicity:** No further relevant information available.
- **Persistence and degradability** No further relevant information available.
- **Behavior in environmental systems:**
- **Bioaccumulative potential** No further relevant information available.
- **Mobility in soil** No further relevant information available.
- **Additional ecological information:**
- **General notes:** Generally not hazardous for water
- **Results of PBT and vPvB assessment**
- **PBT:** Not applicable.
- **vPvB:** Not applicable.
- **Other adverse effects** No further relevant information available.

### 13 Disposal considerations

- **Waste treatment methods**
- **Recommendation:**  
Take to an approved landfill or a waste incineration plant, under conditions approved by the local authority.
- **Uncleaned packagings:**
- **Recommendation:** Disposal must be made according to official regulations.

### 14 Transport information

- **UN-Number**
- **DOT** NA3077
- **ADR, ADN, IMDG, IATA** Void
- **UN proper shipping name**
- **ADR, ADN, IMDG, IATA** Void
- **Transport hazard class(es)**
- **DOT**
- **Class** 9 Miscellaneous dangerous substances and articles
- **Label** 9
- **ADR, ADN, IMDG, IATA**
- **Class** Void
- **Packing group**
- **DOT** III
- **ADR, IMDG, IATA** Void
- **Environmental hazards:**
- **Marine pollutant:** No
- **Special precautions for user** Not applicable.

(Contd. on page 7)

US

## Safety Data Sheet

acc. to OSHA HCS

Printing date 06/23/2015

Version number 5

Reviewed on 05/07/2015

**Trade name: Capricorn**

(Contd. of page 6)

- |  |  |
|--|--|
| <b>· Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code</b> | Not applicable.  |
| <b>· Transport/Additional information:</b>                                       | Product is not classified as a dangerous good for transport (ADR, IMDG, IATA). |
| <b>· UN "Model Regulation":</b>  | -  |

### 15 Regulatory information

- Safety, health and environmental regulations/legislation specific for the substance or mixture

- Sara

- Section 355 (extremely hazardous substances):

None of the ingredients is listed.

- Section 313 (Specific toxic chemical listings):

7440-22-4	silver
7440-50-8	copper

- TSCA (Toxic Substances Control Act):

7440-57-5	Gold
7440-22-4	silver
7440-50-8	copper
7440-06-4	platinum
7440-74-6	indium
7439-88-5	iridium

- Proposition 65

- Chemicals known to cause cancer:

None of the ingredients is listed.

- Chemicals known to cause reproductive toxicity for females:

None of the ingredients is listed.

- Chemicals known to cause reproductive toxicity for males:

None of the ingredients is listed.

- Chemicals known to cause developmental toxicity:

None of the ingredients is listed.

- Carcinogenic categories

- EPA (Environmental Protection Agency)

7440-22-4	silver	
7440-50-8	copper	D

- TLV (Threshold Limit Value established by ACGIH)

None of the ingredients is listed.

- NIOSH-Ca (National Institute for Occupational Safety and Health)

None of the ingredients is listed.

- GHS label elements

The product is classified and labeled according to the Globally Harmonized System (GHS).

(Contd. on page 8)

# Safety Data Sheet

## acc. to OSHA HCS

Printing date 06/23/2015

Version number 5

Reviewed on 05/07/2015

Trade name: **Capricorn**

(Contd. of page 7)

· **Hazard pictograms**

GHS03 GHS05

· **Signal word** *Danger*· **Hazard-determining components of labeling:***gallium*· **Hazard statements***May intensify fire; oxidizer.**Causes severe skin burns and eye damage.*· **Precautionary statements***Keep away from heat.**Take any precaution to avoid mixing with combustibles.**Wear protective gloves / eye protection / face protection.**Wear eye protection / face protection.**Do not breathe dusts or mists.**Keep/Store away from clothing/combustible materials.**Wash thoroughly after handling.**If on skin (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower.**IF INHALED: Remove person to fresh air and keep comfortable for breathing.**If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.**Immediately call a poison center/doctor.**Specific treatment (see on this label).**If swallowed: Rinse mouth. Do NOT induce vomiting.**Wash contaminated clothing before reuse.**In case of fire: Use for extinction: CO<sub>2</sub>, powder or water spray.**Store locked up.**Dispose of contents/container in accordance with local/regional/national/international regulations.*· **National regulations:**· **Other regulations, limitations and prohibitive regulations***The product is a medical device according to the Directive 93/42/EEC.**This product is classified as a medical device under US and Canadian regulations and has been reviewed by the US Food and Drug Administration and Health Canada.*· **Chemical safety assessment:** *A Chemical Safety Assessment has not been carried out.*\* **16 Other information***This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.*· **Date of preparation / last revision** 06/23/2015 / 4· **Abbreviations and acronyms:***ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)**IMDG: International Maritime Code for Dangerous Goods**DOT: US Department of Transportation**IATA: International Air Transport Association**ACGIH: American Conference of Governmental Industrial Hygienists**EINECS: European Inventory of Existing Commercial Chemical Substances**ELINCS: European List of Notified Chemical Substances**CAS: Chemical Abstracts Service (division of the American Chemical Society)**NFPA: National Fire Protection Association (USA)**HMIS: Hazardous Materials Identification System (USA)*

(Contd. on page 9)



**Safety Data Sheet**  
**acc. to OSHA HCS**

Printing date 06/23/2015

Version number 5

Reviewed on 05/07/2015

Trade name: **Capricorn**

*Ox. Sol. 2: Oxidising Solids, Hazard Category 2*  
*Skin Corr. 1C: Skin corrosion/irritation, Hazard Category 1C*  
*Eye Dam. 1: Serious eye damage/eye irritation, Hazard Category 1*  
· **\* Data compared to the previous version altered.**

(Contd. of page 8)

US