

**Safety Data Sheet**  
acc. to OSHA HCS

Printing date 06/23/2015

Version number 5

Reviewed on 05/07/2015


**1 Identification**

- **Product identifier**
- **Trade name:** ***d.SIGN 53***
- **Relevant identified uses of the substance or mixture and uses advised against**  
*No further relevant information available.*
- **Application of the substance / the mixture** *Manufacture of dental prothesis*
- **Details of the supplier of the safety data sheet**
- **Manufacturer/Supplier:**  
*Ivoclar Vivadent Inc.*  
*175 Pineview Drive, Amherst, N.Y. 14228*  
*USA*  
*Tel. +1 800 533 6825*  
*Fax +1 716 691 2285*
  
- Ivoclar Vivadent Inc.*  
*1-6600 Dixie Road*  
*Mississauga, Ontario*  
*L5T 2Y2*  
*Canada*  
*Phone: +1 905 670 8499*  
*Fax: +1 905 670 3102*
- **Information department:** *Quality Assurance / Regulatory Affairs*
- **Emergency telephone number:**  
*24 Hour Emergency Assistance:*  
*Emergency-Call USA - Infotrac: 1-800-535-5053*  
*Emergency-Call Canada - Canutec: 1-613-996-6666*
  
- General SDS Assistance:*  
*US: 1-800-533-6825*  
*Canada: 1-800-263-8182*

**2 Hazard(s) identification**

- **Classification of the substance or mixture**  
*Ox. Sol. 2 H272 May intensify fire; oxidizer.*

---

  - **Label elements**
  - **GHS label elements**  
*The product is classified and labeled according to the Globally Harmonized System (GHS).*
  - **Hazard pictograms**
- 

GHS03
- **Signal word** *Danger*
  - **Hazard statements**  
*May intensify fire; oxidizer.*
  - **Precautionary statements**  
*Keep away from heat.*  
*Take any precaution to avoid mixing with combustibles.*  
*Wear protective gloves / eye protection / face protection.*  
*Keep/Store away from clothing/combustible materials.*

(Contd. on page 2)

# Safety Data Sheet

## acc. to OSHA HCS

Printing date 06/23/2015

Version number 5

Reviewed on 05/07/2015

Trade name: **d.SIGN 53**

(Contd. of page 1)

In case of fire: Use for extinction: CO<sub>2</sub>, powder or water spray.

Dispose of contents/container in accordance with local/regional/national/international regulations.

· **Additional information:**

Alloys do not require a label, providing they do not present a hazard to human health by inhalation, ingestion or contact with skin or to the aquatic environment in the form in which they are placed on the market.

· **Classification system:**

· **NFPA ratings (scale 0 - 4)**



· **HMS-ratings (scale 0 - 4)**

HEALTH	0	Health = 0
FIRE	0	Fire = 0
REACTIVITY	0	Reactivity = 0

· **Other hazards**

· **Results of PBT and vPvB assessment**

· **PBT:** Not applicable.

· **vPvB:** Not applicable.

### 3 Composition/information on ingredients

· **Chemical characterization: Mixtures**

· **Description:** Dental alloy

· **Dangerous components:**

7440-05-3	palladium	50-100%
7440-22-4	silver	25-50%
7440-31-5	tin	2.5-<10%
7440-74-6	indium	1-<2.5%

### 4 First-aid measures

· **Description of first aid measures**

· **General information:** No special measures required.

· **After inhalation:**

Grinding dust: Supply fresh air; consult doctor in case of complaints.

In case of unconsciousness place patient stably in side position for transportation.

· **After skin contact:** After contact with the molten product, cool rapidly with cold water.

· **After eye contact:**

Grinding dust: Rinse opened eye for several minutes under running water. If symptoms persist, consult a doctor.

Mechanical effects only.

· **After swallowing:** Seek medical treatment.

· **Information for doctor:**

· **Most important symptoms and effects, both acute and delayed** No further relevant information available.

· **Indication of any immediate medical attention and special treatment needed**

No further relevant information available.

US

(Contd. on page 3)

# Safety Data Sheet

## acc. to OSHA HCS

Printing date 06/23/2015

Version number 5

Reviewed on 05/07/2015

Trade name: **d.SIGN 53**

(Contd. of page 2)

### 5 Fire-fighting measures

- **Extinguishing media**
- **Suitable extinguishing agents:** *The product is not flammable.*
- **Special hazards arising from the substance or mixture** *No further relevant information available.*
- **Advice for firefighters**
- **Protective equipment:** *Do not inhale explosion gases or combustion gases.*

### 6 Accidental release measures

- **Personal precautions, protective equipment and emergency procedures** *Not required.*
- **Environmental precautions:** *No special measures required.*
- **Methods and material for containment and cleaning up:** *Pick up mechanically.*
- **Reference to other sections**  
*See Section 7 for information on safe handling.*  
*See Section 8 for information on personal protection equipment.*  
*See Section 13 for disposal information.*

### 7 Handling and storage

- **Handling:**
- **Precautions for safe handling**  
*Extractors are required on all machines used for thermal processing or splinter removal processes.*
- **Information about protection against explosions and fires:** *No special measures required.*
- **Conditions for safe storage, including any incompatibilities**
- **Storage:**
- **Requirements to be met by storerooms and receptacles:** *No special requirements.*
- **Information about storage in one common storage facility:** *Not required.*
- **Further information about storage conditions:** *None.*
- **Specific end use(s)** *No further relevant information available.*

### \* 8 Exposure controls/personal protection

- **Additional information about design of technical systems:** *No further data; see item 7.*
- **Control parameters**

· **Components with limit values that require monitoring at the workplace:**

#### 7440-22-4 silver

PEL	Long-term value: 0.01 mg/m <sup>3</sup>
REL	Long-term value: 0.01 mg/m <sup>3</sup>
TLV	Long-term value: 0.1 mg/m <sup>3</sup> metal: dust and fume

#### 7440-31-5 tin

PEL	Long-term value: 2 mg/m <sup>3</sup> metal
REL	Long-term value: 2 mg/m <sup>3</sup>
TLV	Long-term value: 2 mg/m <sup>3</sup> metal

(Contd. on page 4)

US

## Safety Data Sheet

acc. to OSHA HCS

Printing date 06/23/2015

Version number 5

Reviewed on 05/07/2015

**Trade name: d.SIGN 53**

(Contd. of page 3)

### 7440-74-6 indium

REL	Long-term value: 0.1 mg/m <sup>3</sup> as In
TLV	Long-term value: 0.1 mg/m <sup>3</sup> as In

- **Additional information:** The lists that were valid during the creation were used as basis.
- **Exposure controls**
- **Personal protective equipment:**
- **General protective and hygienic measures:** Usual hygienic measures for dental practice.
- **Breathing equipment:** Use respiratory protective device against the effects of fumes/dust/aerosol.
- **Protection of hands:**  
Protective gloves should always be worn during mechanical and thermal processing.  
After use of gloves apply skin-cleaning agents and skin cosmetics.
- **Material of gloves**  
Mechanical processing:  
Leather gloves  
Strong gloves  
Thermal processing:  
Heat protection gloves
- **Penetration time of glove material**  
The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.
- **Eye protection:**  
Always wear safety goggles during mechanical processing (grinding, sawing /cutting, drilling, milling).

## 9 Physical and chemical properties

- **Information on basic physical and chemical properties**
- **General Information**
- **Appearance:**

Form:	Solid
Color:	White
Odor:	Odorless
Odour threshold:	Not determined.
- **pH-value:** Not applicable.
- **Change in condition**

Melting point/Melting range:	1180–1280 °C (2156–2336 °F)
Boiling point/Boiling range:	Undetermined.
- **Flash point:** Not applicable.
- **Flammability (solid, gaseous):** Not determined.
- **Auto igniting:** Product is not selfigniting.
- **Danger of explosion:** Product does not present an explosion hazard.
- **Explosion limits:**

Lower:	Not determined.
Upper:	Not determined.
- **Vapor pressure:** Not applicable.
- **Density at 20 °C (68 °F):** 10.7 g/cm<sup>3</sup> (89.292 lbs/gal)
- **Relative density** Not determined.

(Contd. on page 5)

## Safety Data Sheet

acc. to OSHA HCS

Printing date 06/23/2015

Version number 5

Reviewed on 05/07/2015

Trade name: **d.SIGN 53**

(Contd. of page 4)

· <b>Vapour density</b>	Not applicable.
· <b>Evaporation rate</b>	Not applicable.
· <b>Solubility in / Miscibility with Water:</b>	Insoluble.
· <b>Partition coefficient (n-octanol/water):</b>	Not determined.
· <b>Viscosity:</b>	
<b>Dynamic:</b>	Not applicable.
<b>Kinematic:</b>	Not applicable.
· <b>Solvent content:</b>	
<b>Solids content:</b>	100 %
· <b>Other information</b>	No further relevant information available.

### 10 Stability and reactivity

- **Reactivity** No further relevant information available.
- **Chemical stability** Stable under normal handling and storage conditions.
- **Thermal decomposition / conditions to be avoided:** No decomposition if used according to specifications.
- **Possibility of hazardous reactions** No dangerous reactions known.
- **Conditions to avoid** No further relevant information available.
- **Incompatible materials:** No further relevant information available.
- **Hazardous decomposition products:** None under normal conditions of storage and use.

### 11 Toxicological information

- **Information on toxicological effects**
- **Acute toxicity:**
- **Primary irritant effect:**
- **on the skin:** No irritant effect.
- **on the eye:** No irritating effect.
- **Sensitization:** No sensitizing effects known.
- **Additional toxicological information:**  
Fumes or dusts generated from cutting or grinding operations may cause respiratory irritation.
- **Carcinogenic categories**

#### · **NTP (National Toxicology Program)**

None of the ingredients is listed.

#### · **OSHA-Ca (Occupational Safety & Health Administration)**

None of the ingredients is listed.

### 12 Ecological information

- **Toxicity**
- **Aquatic toxicity:** No further relevant information available.
- **Persistence and degradability** No further relevant information available.
- **Behavior in environmental systems:**
- **Bioaccumulative potential** No further relevant information available.
- **Mobility in soil** No further relevant information available.
- **Additional ecological information:**
- **General notes:** Generally not hazardous for water

(Contd. on page 6)

## Safety Data Sheet acc. to OSHA HCS

Printing date 06/23/2015

Version number 5

Reviewed on 05/07/2015

**Trade name: d.SIGN 53**

(Contd. of page 5)

- **Results of PBT and vPvB assessment**
- **PBT:** Not applicable.
- **vPvB:** Not applicable.
- **Other adverse effects** No further relevant information available.

### 13 Disposal considerations

- **Waste treatment methods**
- **Recommendation:**  
Take to an approved landfill or a waste incineration plant, under conditions approved by the local authority.
- **Uncleaned packagings:**
- **Recommendation:** Disposal must be made according to official regulations.

### 14 Transport information

- |  |  |
|--|--|
| · <b>UN-Number</b>   |  |
| · <b>DOT</b>   | NA3077   |
| · <b>ADR, ADN, IMDG, IATA</b>  | Void   |
| · <b>UN proper shipping name</b>   |  |
| · <b>ADR, ADN, IMDG, IATA</b>  | Void   |
| · <b>Transport hazard class(es)</b>  |  |
| · <b>DOT</b>   |  |
| · <b>Class</b>   | 9 Miscellaneous dangerous substances and articles                              |
| · <b>Label</b>   | 9  |
|  |  |
| · <b>ADR, ADN, IMDG, IATA</b>  |  |
| · <b>Class</b>   | Void   |
| · <b>Packing group</b>   |  |
| · <b>DOT</b>   | III  |
| · <b>ADR, IMDG, IATA</b>   | Void   |
| · <b>Environmental hazards:</b>  |  |
| · <b>Marine pollutant:</b>   | No   |
| · <b>Special precautions for user</b>  | Not applicable.  |
| · <b>Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code</b> | Not applicable.  |
| · <b>Transport/Additional information:</b>                                       | Product is not classified as a dangerous good for transport (ADR, IMDG, IATA). |
| · <b>UN "Model Regulation":</b>  | -  |

### 15 Regulatory information

- **Safety, health and environmental regulations/legislation specific for the substance or mixture**
- **Sara**

· **Section 355 (extremely hazardous substances):**

None of the ingredients is listed.

(Contd. on page 7)

## Safety Data Sheet

acc. to OSHA HCS

Printing date 06/23/2015

Version number 5

Reviewed on 05/07/2015

Trade name: **d.SIGN 53**

(Contd. of page 6)

· **Section 313 (Specific toxic chemical listings):**

7440-22-4	silver
7440-50-8	copper

· **TSCA (Toxic Substances Control Act):**

7440-57-5	Gold
7440-22-4	silver
7440-50-8	copper
7440-06-4	platinum
7440-74-6	indium
7439-88-5	iridium

· **Proposition 65**

· **Chemicals known to cause cancer:**

None of the ingredients is listed.

· **Chemicals known to cause reproductive toxicity for females:**

None of the ingredients is listed.

· **Chemicals known to cause reproductive toxicity for males:**

None of the ingredients is listed.

· **Chemicals known to cause developmental toxicity:**

None of the ingredients is listed.

· **Carcinogenic categories**

· **EPA (Environmental Protection Agency)**

7440-22-4	silver	D
7440-50-8	copper	D

· **TLV (Threshold Limit Value established by ACGIH)**

None of the ingredients is listed.

· **NIOSH-Ca (National Institute for Occupational Safety and Health)**

None of the ingredients is listed.

· **GHS label elements**

The product is classified and labeled according to the Globally Harmonized System (GHS).

· **Hazard pictograms**



GHS03

· **Signal word** *Danger*

· **Hazard statements**

May intensify fire; oxidizer.

· **Precautionary statements**

Keep away from heat.

Take any precaution to avoid mixing with combustibles.

Wear protective gloves / eye protection / face protection.

Keep/Store away from clothing/combustible materials.

In case of fire: Use for extinction: CO<sub>2</sub>, powder or water spray.

Dispose of contents/container in accordance with local/regional/national/international regulations.

(Contd. on page 8)

US

## Safety Data Sheet acc. to OSHA HCS

Printing date 06/23/2015

Version number 5

Reviewed on 05/07/2015

**Trade name: d.SIGN 53**

(Contd. of page 7)

- **National regulations:**
- **Other regulations, limitations and prohibitive regulations**  
The product is a medical device according to the Directive 93/42/EEC.  
This product is classified as a medical device under US and Canadian regulations and has been reviewed by the US Food and Drug Administration and Health Canada.
- **Chemical safety assessment:** A Chemical Safety Assessment has not been carried out.

### \* 16 Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

- **Date of preparation / last revision** 06/23/2015 / 4
- **Abbreviations and acronyms:**  
ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)  
IMDG: International Maritime Code for Dangerous Goods  
DOT: US Department of Transportation  
IATA: International Air Transport Association  
ACGIH: American Conference of Governmental Industrial Hygienists  
EINECS: European Inventory of Existing Commercial Chemical Substances  
ELINCS: European List of Notified Chemical Substances  
CAS: Chemical Abstracts Service (division of the American Chemical Society)  
NFPA: National Fire Protection Association (USA)  
HMIS: Hazardous Materials Identification System (USA)  
Ox. Sol. 2: Oxidising Solids, Hazard Category 2
- **\* Data compared to the previous version altered.**

US