

**Safety Data Sheet**  
acc. to OSHA HCS

Printing date 06/23/2015

Version number 7

Reviewed on 05/07/2015


**1 Identification**

- **Product identifier**
- **Trade name: d.SIGN 96**
- **Relevant identified uses of the substance or mixture and uses advised against**  
No further relevant information available.
- **Application of the substance / the mixture** Manufacture of dental prothesis
- **Details of the supplier of the safety data sheet**
- **Manufacturer/Supplier:**  
Ivoclar Vivadent Inc.  
175 Pineview Drive, Amherst, N.Y. 14228  
USA  
Tel. +1 800 533 6825  
Fax +1 716 691 2285
  
- Ivoclar Vivadent Inc.  
1-6600 Dixie Road  
Mississauga, Ontario  
L5T 2Y2  
Canada  
Phone: +1 905 670 8499  
Fax: +1 905 670 3102
- **Information department:** Quality Assurance / Regulatory Affairs
- **Emergency telephone number:**  
24 Hour Emergency Assistance:  
Emergency-Call USA - Infotrac: 1-800-535-5053  
Emergency-Call Canada - Canutec: 1-613-996-6666
  
- General SDS Assistance:  
US: 1-800-533-6825  
Canada: 1-800-263-8182

\* **2 Hazard(s) identification**

- **Classification of the substance or mixture**  
Ox. Sol. 2 H272 May intensify fire; oxidizer.

---

  - **Label elements**
  - **GHS label elements**  
The product is classified and labeled according to the Globally Harmonized System (GHS).
  - **Hazard pictograms**
- 

GHS03
- **Signal word** Danger
  - **Hazard statements**  
May intensify fire; oxidizer.
  - **Precautionary statements**  
Keep away from heat.  
Take any precaution to avoid mixing with combustibles.  
Wear protective gloves / eye protection / face protection.  
Keep/Store away from clothing/combustible materials.

(Contd. on page 2)

# Safety Data Sheet

## acc. to OSHA HCS

Printing date 06/23/2015

Version number 7

Reviewed on 05/07/2015

**Trade name: d.SIGN 96**

(Contd. of page 1)

*In case of fire: Use for extinction: CO<sub>2</sub>, powder or water spray.*

*Dispose of contents/container in accordance with local/regional/national/international regulations.*

· **Additional information:**

*Alloys do not require a label, providing they do not present a hazard to human health by inhalation, ingestion or contact with skin or to the aquatic environment in the form in which they are placed on the market.*

· **Information pertaining to particular dangers for man and environment:**

*Alloys do not require a label, providing they do not present a hazard to human health by inhalation, ingestion or contact with skin or to the aquatic environment in the form in which they are placed on the market.*

· **Classification system:**

· **NFPA ratings (scale 0 - 4)**



· **HMIS-ratings (scale 0 - 4)**



· **Other hazards**

· **Results of PBT and vPvB assessment**

· **PBT:** Not applicable.

· **vPvB:** Not applicable.

### 3 Composition/information on ingredients

· **Chemical characterization: Mixtures**

· **Description:** Mixture of the substances listed below with nonhazardous additions.

· **Dangerous components:**

|           |           |          |
|-----------|-----------|----------|
| 7440-22-4 | silver    | 2.5-<10% |
| 7440-06-4 | platinum  | 2.5-<10% |
| 7440-05-3 | palladium | 2.5-<10% |
| 7440-74-6 | indium    | 1-<2.5%  |

### 4 First-aid measures

· **Description of first aid measures**

· **After inhalation:** Grinding dust: Supply fresh air; consult doctor in case of complaints.

· **After skin contact:** After contact with the molten product, cool rapidly with cold water.

· **After eye contact:**

*Grinding dust: Rinse opened eye for several minutes under running water. If symptoms persist, consult a doctor.*

*Mechanical effects only.*

· **After swallowing:** Seek medical treatment.

· **Information for doctor:**

· **Most important symptoms and effects, both acute and delayed** No further relevant information available.

· **Indication of any immediate medical attention and special treatment needed**

*No further relevant information available.*

US

(Contd. on page 3)

## Safety Data Sheet acc. to OSHA HCS

Printing date 06/23/2015

Version number 7

Reviewed on 05/07/2015

Trade name: **d.SIGN 96**

(Contd. of page 2)

### 5 Fire-fighting measures

- **Extinguishing media**
- **Suitable extinguishing agents:** *The product is not flammable.*
- **Special hazards arising from the substance or mixture** *No further relevant information available.*
- **Advice for firefighters**
- **Protective equipment:** *Do not inhale explosion gases or combustion gases.*

### 6 Accidental release measures

- **Personal precautions, protective equipment and emergency procedures**  
*Wear protective equipment. Keep unprotected persons away.*
- **Environmental precautions:** *No special measures required.*
- **Methods and material for containment and cleaning up:** *Pick up mechanically.*
- **Reference to other sections**  
*See Section 7 for information on safe handling.*  
*See Section 8 for information on personal protection equipment.*  
*See Section 13 for disposal information.*

### 7 Handling and storage

- **Handling:**
- **Precautions for safe handling**  
*Extractors are required on all machines used for thermal processing or splinter removal processes.*  
*Only adequately trained personnel should handle this product.*  
*For use in dentistry only.*
- **Information about protection against explosions and fires:** *Keep ignition sources away - Do not smoke.*
- **Conditions for safe storage, including any incompatibilities**
- **Storage:**
- **Requirements to be met by storerooms and receptacles:** *Store only in the original receptacle.*
- **Information about storage in one common storage facility:** *Not required.*
- **Further information about storage conditions:** *None.*
- **Specific end use(s)** *No further relevant information available.*

### \* 8 Exposure controls/personal protection

- **Additional information about design of technical systems:** *No further data; see item 7.*
- **Control parameters**

· **Components with limit values that require monitoring at the workplace:**

#### 7440-22-4 silver

PEL Long-term value: 0.01 mg/m<sup>3</sup>

REL Long-term value: 0.01 mg/m<sup>3</sup>

TLV Long-term value: 0.1 mg/m<sup>3</sup>  
metal: dust and fume

#### 7440-06-4 platinum

PEL Long-term value: 0.002 mg/m<sup>3</sup>  
soluble salts, as Pt

(Contd. on page 4)

US

## Safety Data Sheet

### acc. to OSHA HCS

Printing date 06/23/2015

Version number 7

Reviewed on 05/07/2015

**Trade name: d.SIGN 96**

(Contd. of page 3)

REL Long-term value: 1\* 0.002\*\* mg/m<sup>3</sup>  
\*Metal \*\*soluble salts, as Pt

TLV Long-term value: 1\* 0.002\*\* mg/m<sup>3</sup>  
\*metal \*\*soluble salts, as Pt

#### 7440-74-6 indium

REL Long-term value: 0.1 mg/m<sup>3</sup>  
as In

TLV Long-term value: 0.1 mg/m<sup>3</sup>  
as In

· **Additional information:** The lists that were valid during the creation were used as basis.

· **Exposure controls**

· **Personal protective equipment:**

· **General protective and hygienic measures:**

Wash hands before breaks and at the end of work.

Usual hygienic measures for dental practice.

· **Breathing equipment:** Use respiratory protective device against the effects of fumes/dust/aerosol.

· **Protection of hands:**

Protective gloves should always be worn during mechanical and thermal processing.

After use of gloves apply skin-cleaning agents and skin cosmetics.

· **Material of gloves**

Mechanical processing:

Leather gloves

Strong gloves

Thermal processing:

Heat protection gloves

· **Penetration time of glove material**

The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.

· **Eye protection:**

Always wear safety goggles during mechanical processing (grinding, sawing /cutting, drilling, milling).

## 9 Physical and chemical properties

· **Information on basic physical and chemical properties**

· **General Information**

· **Appearance:**

Form: Solid

Color: Yellow

· **Odor:** Odorless

· **Odour threshold:** Not determined.

· **pH-value:** Not applicable.

· **Change in condition**

Melting point/Melting range: 1050–1170 °C (1922–2138 °F)

Boiling point/Boiling range: Undetermined.

· **Flash point:** Not applicable.

· **Flammability (solid, gaseous):** Contact with combustible material may cause fire.

· **Auto igniting:** Product is not selfigniting.

· **Danger of explosion:** Not determined.

(Contd. on page 5)

## Safety Data Sheet

acc. to OSHA HCS

Printing date 06/23/2015

Version number 7

Reviewed on 05/07/2015

Trade name: **d.SIGN 96**

(Contd. of page 4)

|   |  |
|---|--|
| · <b>Explosion limits:</b>                        |  |
| <b>Lower:</b>                                     | Not determined.                            |
| <b>Upper:</b>                                     | Not determined.                            |
| · <b>Vapor pressure:</b>                          | Not applicable.                            |
| · <b>Density at 20 °C (68 °F):</b>                | 16.7 g/cm <sup>3</sup> (139.362 lbs/gal)   |
| · <b>Relative density</b>                         | Not determined.                            |
| · <b>Vapour density</b>                           | Not applicable.                            |
| · <b>Evaporation rate</b>                         | Not applicable.                            |
| · <b>Solubility in / Miscibility with Water:</b>  | Insoluble.                                 |
| · <b>Partition coefficient (n-octanol/water):</b> | Not determined.                            |
| · <b>Viscosity:</b>                               |  |
| <b>Dynamic:</b>                                   | Not applicable.                            |
| <b>Kinematic:</b>                                 | Not applicable.                            |
| · <b>Solvent content:</b>                         |  |
| <b>Organic solvents:</b>                          | 0.0 %                                      |
| <b>Solids content:</b>                            | 100.0 %                                    |
| · <b>Other information</b>                        | No further relevant information available. |

### 10 Stability and reactivity

- **Reactivity** No further relevant information available.
- **Chemical stability** Stable under normal handling and storage conditions.
- **Thermal decomposition / conditions to be avoided:** No decomposition if used according to specifications.
- **Possibility of hazardous reactions** No dangerous reactions known.
- **Conditions to avoid** No further relevant information available.
- **Incompatible materials:** No further relevant information available.
- **Hazardous decomposition products:** None under normal conditions of storage and use.

### 11 Toxicological information

- **Information on toxicological effects**
- **Acute toxicity:**
- **Primary irritant effect:**
- **on the skin:** No irritant effect.
- **on the eye:** No irritating effect.
- **Sensitization:** No sensitizing effects known.
- **Additional toxicological information:**  
Fumes or dusts generated from cutting or grinding operations may cause respiratory irritation.
- **Carcinogenic categories**

#### · **NTP (National Toxicology Program)**

None of the ingredients is listed.

#### · **OSHA-Ca (Occupational Safety & Health Administration)**

None of the ingredients is listed.

US

(Contd. on page 6)

## Safety Data Sheet

acc. to OSHA HCS

Printing date 06/23/2015

Version number 7

Reviewed on 05/07/2015

Trade name: **d.SIGN 96**

(Contd. of page 5)

### 12 Ecological information

- **Toxicity**
- **Aquatic toxicity:** No further relevant information available.
- **Persistence and degradability** No further relevant information available.
- **Behavior in environmental systems:**
- **Bioaccumulative potential** No further relevant information available.
- **Mobility in soil** No further relevant information available.
- **Additional ecological information:**
- **General notes:** Generally not hazardous for water
- **Results of PBT and vPvB assessment**
- **PBT:** Not applicable.
- **vPvB:** Not applicable.
- **Other adverse effects** No further relevant information available.

### 13 Disposal considerations

- **Waste treatment methods**
- **Recommendation:**  
Must not be disposed of together with household garbage. Do not allow product to reach sewage system.  
Take to an approved landfill or a waste incineration plant, under conditions approved by the local authority.
- **Uncleaned packagings:**
- **Recommendation:** Disposal must be made according to official regulations.

### 14 Transport information

- |   |  |
|---|--|
| · <b>UN-Number</b><br>· <b>DOT, ADR, ADN, IMDG, IATA</b>                                    | Void   |
| · <b>UN proper shipping name</b><br>· <b>DOT, ADR, ADN, IMDG, IATA</b>                      | Void   |
| · <b>Transport hazard class(es)</b><br>· <b>DOT, ADR, ADN, IMDG, IATA</b><br>· <b>Class</b> | Void   |
| · <b>Packing group</b><br>· <b>DOT, ADR, IMDG, IATA</b>                                     | Void   |
| · <b>Environmental hazards:</b><br>· <b>Marine pollutant:</b>                               | No   |
| · <b>Special precautions for user</b>   | Not applicable.  |
| · <b>Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code</b>            | Not applicable.  |
| · <b>Transport/Additional information:</b>  | Product is not classified as a dangerous good for transport (ADR, IMDG, IATA). |
| · <b>UN "Model Regulation":</b>   | -  |

US  
(Contd. on page 7)

# Safety Data Sheet

## acc. to OSHA HCS

Printing date 06/23/2015

Version number 7

Reviewed on 05/07/2015

Trade name: **d.SIGN 96**

(Contd. of page 6)

### 15 Regulatory information

· **Safety, health and environmental regulations/legislation specific for the substance or mixture**

· **Sara**

· **Section 355 (extremely hazardous substances):**

None of the ingredients is listed.

· **Section 313 (Specific toxic chemical listings):**

|           |        |
|-----------|--------|
| 7440-22-4 | silver |
|-----------|--------|

|           |           |
|-----------|-----------|
| 7439-96-5 | manganese |
|-----------|-----------|

· **TSCA (Toxic Substances Control Act):**

All ingredients are listed.

· **Proposition 65**

· **Chemicals known to cause cancer:**

None of the ingredients is listed.

· **Chemicals known to cause reproductive toxicity for females:**

None of the ingredients is listed.

· **Chemicals known to cause reproductive toxicity for males:**

None of the ingredients is listed.

· **Chemicals known to cause developmental toxicity:**

None of the ingredients is listed.

· **Carcinogenic categories**

· **EPA (Environmental Protection Agency)**

|           |        |   |
|-----------|--------|---|
| 7440-22-4 | silver | D |
|-----------|--------|---|

|           |           |   |
|-----------|-----------|---|
| 7439-96-5 | manganese | D |
|-----------|-----------|---|

· **TLV (Threshold Limit Value established by ACGIH)**

None of the ingredients is listed.

· **NIOSH-Ca (National Institute for Occupational Safety and Health)**

None of the ingredients is listed.

· **GHS label elements**

The product is classified and labeled according to the Globally Harmonized System (GHS).

· **Hazard pictograms**



GHS03

· **Signal word** *Danger*

· **Hazard statements**

May intensify fire; oxidizer.

· **Precautionary statements**

Keep away from heat.

Take any precaution to avoid mixing with combustibles.

Wear protective gloves / eye protection / face protection.

Keep/Store away from clothing/combustible materials.

In case of fire: Use for extinction: CO<sub>2</sub>, powder or water spray.

Dispose of contents/container in accordance with local/regional/national/international regulations.

(Contd. on page 8)

## Safety Data Sheet acc. to OSHA HCS

Printing date 06/23/2015

Version number 7

Reviewed on 05/07/2015

**Trade name: d.SIGN 96**

(Contd. of page 7)

- **National regulations:**
- **Other regulations, limitations and prohibitive regulations**  
The product is a medical device according to the Directive 93/42/EEC.  
This product is classified as a medical device under US and Canadian regulations and has been reviewed by the US Food and Drug Administration and Health Canada.
- **Chemical safety assessment:** A Chemical Safety Assessment has not been carried out.

### \* 16 Other information

*This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.*

- **Date of preparation / last revision** 06/23/2015 / 6
- **Abbreviations and acronyms:**  
ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)  
IMDG: International Maritime Code for Dangerous Goods  
DOT: US Department of Transportation  
IATA: International Air Transport Association  
ACGIH: American Conference of Governmental Industrial Hygienists  
EINECS: European Inventory of Existing Commercial Chemical Substances  
ELINCS: European List of Notified Chemical Substances  
CAS: Chemical Abstracts Service (division of the American Chemical Society)  
NFPA: National Fire Protection Association (USA)  
HMIS: Hazardous Materials Identification System (USA)  
Ox. Sol. 2: Oxidising Solids, Hazard Category 2
- **\* Data compared to the previous version altered.**

US