



**Skyce**<sup>®</sup>  
Tooth jewels

**A touch  
of glamour**

**ivoclar**  
**vivadent**<sup>®</sup>  
passion vision innovation

# Create sparkling accents

The dazzling Skyce® crystals add a flash of extra sparkle to a smile. The tooth jewels come with a coordinated adhesive system so that you can provide the high-quality professional results that your patients expect.

**Benefits** for your customers:

**Preservation of tooth structure** due to non-invasive preparation and bonding procedure

**Coordinated product system** with Total Etch and Heliobond to ensure a reliable bond and a gentle residue-free removal

**Available in two sizes:**  
ø 1.9 mm and ø 2.4 mm

**Magical light effects**  
due to premium quality crystal glass



Etch dental enamel with phosphoric acid



Apply Heliobond adhesive



Place Skyce® onto the bonding site



Light cure from all sides

Source:  
Dr A. Peschke, Ivoclar Vivadent AG,  
Schaan, Liechtenstein

Art. no.	Description
559320AN	Skyce® Intro Package
559322AN	Skyce® Refill 5 Crystal, large DE
559321AN	Skyce® Refill 5 Crystal, small DE



Production and Distribution  
**Ivoclar Vivadent AG**  
Benderstr. 2  
9494 Schaan  
Liechtenstein  
Tel. +423 235 35 35  
Fax +423 235 33 60  
[www.ivoclarvivadent.com](http://www.ivoclarvivadent.com)

**ivoclar**  
**vivadent**<sup>®</sup>  
passion vision innovation