

**Safety Data Sheet**  
according to WHS Regulations

Printing date 06.12.2019

Version number 1

Revision: 06.12.2019

Hazardous according to criteria of Australian Safety and Compensation Council.

## 1 Identification

- **Product identifier**
- **Trade name:** ***IvoBase CAD Bond Polymer***
- **Relevant identified uses of the substance or mixture and uses advised against**  
No further relevant information available.
- **Application of the substance / the mixture** Auxiliary for manufacture of dental prothesis
- **Details of the supplier of the safety data sheet**
- **Manufacturer/Supplier:**  
Ivoclar Vivadent AG  
Bendererstrasse 2  
9494 Schaan  
PRINCIPALITY OF LIECHTENSTEIN  
Tel: +423 235 35 35 / Fax: +423 235 33 60
- **Importer:**  
Ivoclar Vivadent Pty. Ltd.  
1- 5 Overseas Drive  
Noble Park North VIC 3174  
Tel: + 61 3 9795 9599 / Fax: + 61 3 9795 9645
- **Further information obtainable from:**  
Regulatory Affairs  
sds@ivoclarvivadent.com
- **Emergency telephone number:** 131 126 (Poisons Information Centre - 24 hours / 7 days)

## 2 Hazard(s) Identification

- **Classification of the substance or mixture**  
Skin Sens. 1 H317 May cause an allergic skin reaction.
- **Label elements**
- **GHS label elements**  
The product is classified and labelled according to the Globally Harmonised System (GHS).
- **Hazard pictograms**



GHS07

- **Signal word** Warning
- **Hazard-determining components of labelling:**  
dibenzoyl peroxide
- **Hazard statements**  
May cause an allergic skin reaction.
- **Precautionary statements**  
Avoid breathing dust/fume/gas/mist/vapours/spray.  
Do not get in eyes, on skin, or on clothing.  
Wear protective gloves/protective clothing/eye protection/face protection.  
Dispose of contents/container in accordance with local/regional/national/international regulations.

(Contd. on page 2)

## Safety Data Sheet according to WHS Regulations

Printing date 06.12.2019

Version number 1

Revision: 06.12.2019

**Trade name: IvoBase CAD Bond Polymer**

(Contd. of page 1)

- **Other hazards**
  - Dust generation
  - Particular danger of slipping on leaked/spilled product.
- **Results of PBT and vPvB assessment**
- **PBT:** Not applicable.
- **vPvB:** Not applicable.

### 3 Composition and Information on Ingredients

- **Chemical characterisation: Mixtures**
- **Description:** > 95% polymethylmethacrylate

#### · **Dangerous components:**

CAS: 94-36-0	dibenzoyl peroxide	Org. Perox. B, H241; Eye Irrit. 2, H319; Skin Sens. 1, H317	<2.5%
--------------	--------------------	---	-------

- **Additional information:** For the wording of the listed hazard phrases refer to section 16.

### 4 First Aid Measures

- **Description of first aid measures**
- **General information:** No special measures required.
- **After inhalation:** Supply fresh air and to be sure call for a doctor.
- **After skin contact:**
  - Rinse with water.
  - If skin irritation continues, consult a doctor.
- **After eye contact:** Rinse opened eye for several minutes under running water. Then consult a doctor.
- **After swallowing:**
  - Rinse out mouth and then drink plenty of water.
  - If symptoms persist consult doctor.
- **Information for doctor:**
- **Most important symptoms and effects, both acute and delayed** No further relevant information available.
- **Indication of any immediate medical attention and special treatment needed**
  - No further relevant information available.

### 5 Fire Fighting Measures

- **Extinguishing media**
- **Suitable extinguishing agents:** Use fire extinguishing methods suitable to surrounding conditions.
- **Special hazards arising from the substance or mixture** No further relevant information available.
- **Advice for firefighters**
- **Protective equipment:** No special measures required.

### 6 Accidental Release Measures

- **Personal precautions, protective equipment and emergency procedures** Not required.
- **Environmental precautions:** No special measures required.
- **Methods and material for containment and cleaning up:** Pick up mechanically.
- **Reference to other sections**
  - See Section 7 for information on safe handling.
  - See Section 8 for information on personal protection equipment.
  - See Section 13 for disposal information.

AU

(Contd. on page 3)

## Safety Data Sheet according to WHS Regulations

Printing date 06.12.2019

Version number 1

Revision: 06.12.2019

**Trade name: IvoBase CAD Bond Polymer**

(Contd. of page 2)

### 7 Handling and Storage

- **Handling:**
- **Precautions for safe handling**  
Only adequately trained personnel should handle this product.  
For use in dentistry only.
- **Information about fire - and explosion protection:** No special measures required.
- **Conditions for safe storage, including any incompatibilities**
- **Storage:**
- **Requirements to be met by storerooms and receptacles:** Store only in the original receptacle.
- **Information about storage in one common storage facility:** Not required.
- **Further information about storage conditions:** Protect from heat and direct sunlight.
- **Specific end use(s)** No further relevant information available.

### 8 Exposure controls and personal protection

- **Additional information about design of technical facilities:** No further data; see item 7.
- **Control parameters**

· **Ingredients with limit values that require monitoring at the workplace:**

**CAS: 94-36-0 dibenzoyl peroxide**

NES	Long-term value: 5 mg/m <sup>3</sup> Sen
WES	Long-term value: 5 mg/m <sup>3</sup> Sen

- **Additional information:** The lists valid during the making were used as basis.
- **Exposure controls**
- **Personal protective equipment:**
- **General protective and hygienic measures:**  
Usual hygienic measures for dental practice and dental laboratories.  
Keep away from foodstuffs, beverages and feed.  
Wash hands before breaks and at the end of work.  
Avoid contact with the skin.
- **Respiratory protection:** Use respiratory protective device against the effects of dust.
- **Protection of hands:**



Protective gloves

- **Material of gloves**  
Butyl rubber, BR  
Natural rubber, NR  
Nitrile rubber, NBR  
Fluorocarbon rubber (Viton)  
Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation
- **Penetration time of glove material**  
The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.

(Contd. on page 4)

## Safety Data Sheet

### according to WHS Regulations

Printing date 06.12.2019

Version number 1

Revision: 06.12.2019

**Trade name: IvoBase CAD Bond Polymer**

(Contd. of page 3)

· **Eye protection:** Safety glasses

### 9 Physical and Chemical Properties

· **Information on basic physical and chemical properties**

· **General Information**

· **Appearance:**

· **Form:** Powder  
 · **Colour:** According to product specification  
 · **Odour:** Characteristic  
 · **Odour threshold:** Not determined.

· **pH-value:** Not applicable.

· **Change in condition**

· **Melting point/freezing point:** 150 °C  
 · **Initial boiling point and boiling range:** Undetermined.

· **Flash point:** 250 °C

· **Flammability (solid, gas):** Not determined.

· **Auto-ignition temperature:** Product is not selfigniting.

· **Explosive properties:** Product does not present an explosion hazard.

· **Vapour pressure:** Not applicable.

· **Density at 20 °C:** 1.2 g/cm<sup>3</sup>  
 · **Relative density:** Not determined.  
 · **Vapour density:** Not applicable.  
 · **Evaporation rate:** Not applicable.

· **Solubility in / Miscibility with water:** Nearly insoluble.

· **Partition coefficient: n-octanol/water:** Not determined.

· **Viscosity:**

· **Dynamic:** Not applicable.  
 · **Kinematic:** Not applicable.

· **Other information:** No further relevant information available.

### 10 Stability and Reactivity

- **Reactivity:** No further relevant information available.
- **Chemical stability:** Stable under normal handling and storage conditions.
- **Thermal decomposition / conditions to be avoided:** No decomposition if used according to specifications.
- **Possibility of hazardous reactions:** No dangerous reactions known.
- **Conditions to avoid:** No further relevant information available.
- **Incompatible materials:** No further relevant information available.
- **Hazardous decomposition products:** None under normal conditions of storage and use.

-AU-

(Contd. on page 5)

## Safety Data Sheet according to WHS Regulations

Printing date 06.12.2019

Version number 1

Revision: 06.12.2019

**Trade name: IvoBase CAD Bond Polymer**

(Contd. of page 4)

### 11 Toxicological Information

- **Information on toxicological effects**
- **Acute toxicity**
- **Skin corrosion/irritation** No irritant effect.
- **Serious eye damage/irritation** No irritating effect.
- **Respiratory or skin sensitisation** Sensitisation possible through skin contact.
- **Additional toxicological information:** No further relevant information available.

### 12 Ecological Information

- **Toxicity**
- **Aquatic toxicity:** No further relevant information available.
- **Persistence and degradability** No further relevant information available.
- **Behaviour in environmental systems:**
- **Bioaccumulative potential** No further relevant information available.
- **Mobility in soil** No further relevant information available.
- **Additional ecological information:**
- **General notes:** Not hazardous for water.
- **Results of PBT and vPvB assessment**
- **PBT:** Not applicable.
- **vPvB:** Not applicable.
- **Other adverse effects** No further relevant information available.

### 13 Disposal considerations

- **Waste treatment methods**
- **Recommendation**  
Take to an approved landfill or a waste incineration plant, under conditions approved by the local authority.
- **Uncleaned packaging:**
- **Recommendation:** Disposal must be made according to official regulations.

### 14 Transport information

- |                                       |                 |
|---------------------------------------|-----------------|
| · <b>UN-Number</b>                    |                 |
| · <b>ADG, ADN, IMDG, IATA</b>         | Void            |
| · <b>UN proper shipping name</b>      |                 |
| · <b>ADG, ADN, IMDG, IATA</b>         | Void            |
| · <b>Transport hazard class(es)</b>   |                 |
| · <b>ADG, ADN, IMDG, IATA</b>         |                 |
| · <b>Class</b>                        | Void            |
| · <b>Packing group</b>                |                 |
| · <b>ADG, IMDG, IATA</b>              | Void            |
| · <b>Environmental hazards:</b>       |                 |
| · <b>Marine pollutant:</b>            | No              |
| · <b>Special precautions for user</b> | Not applicable. |

(Contd. on page 6)

-AU-

## Safety Data Sheet

according to WHS Regulations

Printing date 06.12.2019

Version number 1

Revision: 06.12.2019

**Trade name: IvoBase CAD Bond Polymer**

(Contd. of page 5)

- |   |  |
|---|--|
| <b>· Transport in bulk according to Annex II of Marpol and the IBC Code</b> | Not applicable.  |
| <b>· Transport/Additional information:</b>                                  | Product is not classified as a dangerous good for transport (ADR, IMDG, IATA). |
| <b>· UN "Model Regulation":</b>   | Void   |

### 15 Regulatory information

- **Safety, health and environmental regulations/legislation specific for the substance or mixture**

- **Australian Inventory of Chemical Substances**

All ingredients are listed.

- **Standard for the Uniform Scheduling of Medicines and Poisons**

CAS: 94-36-0 dibenzoyl peroxide

S2, S4, S5

- **GHS label elements**

The product is classified and labelled according to the Globally Harmonised System (GHS).

- **Hazard pictograms**



GHS07

- **Signal word** Warning

- **Hazard-determining components of labelling:**

dibenzoyl peroxide

- **Hazard statements**

May cause an allergic skin reaction.

- **Precautionary statements**

Avoid breathing dust/fume/gas/mist/vapours/spray.

Do not get in eyes, on skin, or on clothing.

Wear protective gloves/protective clothing/eye protection/face protection.

Dispose of contents/container in accordance with local/regional/national/international regulations.

- **Directive 2012/18/EU**

- **Named dangerous substances - ANNEX I** None of the ingredients is listed.

- **National regulations:**

- **Other regulations, limitations and prohibitive regulations**

The product is a medical device according to the Directive 93/42/EEC.

- **Chemical safety assessment:** A Chemical Safety Assessment has not been carried out.

### 16 Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

- **Relevant phrases**

H241 Heating may cause a fire or explosion.

H317 May cause an allergic skin reaction.

H319 Causes serious eye irritation.

(Contd. on page 7)

AU

## **Safety Data Sheet** according to WHS Regulations

Printing date 06.12.2019

Version number 1

Revision: 06.12.2019

**Trade name: IvoBase CAD Bond Polymer**

(Contd. of page 6)

**Abbreviations and acronyms:**

*ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)*

*IMDG: International Maritime Code for Dangerous Goods*

*IATA: International Air Transport Association*

*EINECS: European Inventory of Existing Commercial Chemical Substances*

*ELINCS: European List of Notified Chemical Substances*

*CAS: Chemical Abstracts Service (division of the American Chemical Society)*

*PBT: Persistent, Bioaccumulative and Toxic*

*vPvB: very Persistent and very Bioaccumulative*

*Org. Perox. B: Organic peroxides – Type B*

*Eye Irrit. 2: Serious eye damage/eye irritation – Category 2*

*Skin Sens. 1: Skin sensitisation – Category 1*

-AU-