

IPS e.max® ZirCAD MT Multi

Combination table / Kombinationstabelle

All ceramic,
all you need.

Staining technique (monolithic restorations)

Maltechnik (monolithische Restaurationen)

Tooth shade Zahnfarbe	IPS e.max ZirCAD MT Multi	IPS Ivocolor Shade Dentin	IPS Ivocolor Shade Incisal
BL1	BL1	0	1/2/3
BL2	BL1	0	
BL3	BL3	0	
BL4	BL3	0	
A1	A1	1	
A2	A2	2	
A3	A3	2	
A3.5	A3.5	2	
A4	A3.5	2	
B1	B1	4	
B2	B2	4	
B3	B2	5	
B4	A3	5	
C1	C1	6	
C2	C2	6	
C3	C3	6	
C4	C3	6	
D2	D2	7	
D3	D3	7	
D4	C2	6	



see instructions

For further details on the processing of IPS e.max ZirCAD please refer to the instruction for use IPS e.max ZirCAD Labside.

Für weitere Verarbeitungsdetails von IPS e.max ZirCAD beachten Sie bitte die IPS e.max ZirCAD Labside Gebrauchsinformationen

CE 0123



Rx ONLY
For dental use only

Manufacturer:
Ivoclar Vivadent AG, 9494 Schaan/Liechtenstein
www.ivoclarvivadent.com

Date information prepared:
2020-01-31, Rev. 0

ivoclar
vivadent®
passion vision innovation