

# Programat® P300/G2

## Software Revision Log

Version	Changes from previous version	Status
<b>V8.1</b> 25.08.2017	<b>Software optimized</b> <ul style="list-style-type: none"> <li>New program structure</li> <li>New programs (e.g. IPS e.max Ceram)</li> </ul>	<b>Currently delivered version</b>
<b>V8.0</b> 01.04.2016	<b>Software optimized</b> <ul style="list-style-type: none"> <li>New programs for IPS Style and IPS Ivocolor</li> </ul>	<b>No longer available!</b>
<b>V7.1</b> 01.10.2014	<b>Software optimized</b> <ul style="list-style-type: none"> <li>New programs for IPS InLine powder opaquer</li> <li>Various software improvements</li> </ul>	<b>No longer available!</b>
<b>V7.0</b> 05.05.2014	<b>Software optimized</b> <ul style="list-style-type: none"> <li>New service interval indicator</li> <li>Various software improvements</li> <li>New, special speed crystallization program for IPS e.max CAD (P85): The new program replaces the existing IPS e.max CAD Speed Crystallization. The program may only be used in conjunction with the IPS e.max CAD Speed Crystallization Tray (article code 647633AN). Please observe the processing information for the IPS e.max CAD Speed Crystallization.</li> </ul>	<b>No longer available!</b>
<b>V6.0</b> 04.04.2011	<b>Software optimized</b> <ul style="list-style-type: none"> <li>New Ivoclar Vivadent Programs for IPS InLine One (Programs 121 – 126)</li> </ul>	<b>No longer available!</b>
<b>V5.2</b> 28.01.2011	<b>Software optimized</b> <ul style="list-style-type: none"> <li>Update function has been optimized</li> </ul>	<b>No longer available!</b>
<b>V5.1</b> 10.09.2010	<b>Software optimized</b> <ul style="list-style-type: none"> <li>Integration of the Ivoclar Vivadent Programs for the "IPS e.max CAD-on technique" (P83 and P84)</li> </ul>	<b>No longer available!</b>
<b>V5.0</b> 02.08.2010	<b>Basic Version of the Programat P300/G2</b> Compared to the predecessor model Programat P300 following functions got integrated: <ul style="list-style-type: none"> <li>Power Saving Technology</li> <li>TSP: Thermo Shock Protection</li> <li>New program-structure: 120 Ivoclar Vivadent and 300 individual programs</li> </ul>	<b>No longer available!</b>

Version 25.08.2017