

Telio®

Temporary solution out of one hand

ivoclar
vivadent®
clinical

Telio CAD



Telio CS C&B



Activation
Aktivierung



Bonding
Haftvermittler



Individualization
Individualisierung



Modifications, Additions
Modifikationen, Ergänzungen



Polishing
Politur



Desensitization
Desensibilisierung




Cementation
Befestigung



Telio CS Inlay/Onlay



Modifications, Additions
Modifikationen, Ergänzungen



Relining
Unterfütterung



Individualization
Individualisierung



CS C&B

CAD

CS C&B

CAD

CS C&B

CAD

Telio Activator
Activation / Aktivierung

X

X

X

Heliobond
Bonding / Haftvermittler

(X)*

X

X

X

X

Telio Add-On Flow
Modifications, Additions
Modifikationen, Ergänzungen

X

X

Telio CS C&B

X

X

Telio Stains (blue, light-yellow, white, orange)
Individualization / Individualisierung

X

X

OptraPol NG
Polishing / Politur

X

X

X

X

X

X

Telio CS Desensitizer (optional)
Desensitization / Desensibilisierung

(X)

(X)

(X)

(X)

(X)

(X)

Telio CS Link
Cementation / Befestigung

X

X

X

X

X

X