

Telio®

Temporary solution out of one hand

ivoclar
vivadent®
technical

Telio CAD



Telio Lab



Activation Aktivierung



Bonding Haftvermittler



Individualization Individualisierung



Relining / Cut-Back Unterfütterung / Cut-Back



Modifications / Emergence profile Modifikationen / Emergenzprofil



Cementation Befestigung



Prevention of inhibition layer Verhinderung der Inhibitionsschicht




Individualization / Cut-Back
Individualisierung / Cut-Back

Relining
Unterfütterung

Modifications
Modifikationen

Emergence profile
Emergenzprofil

Occlusion build up
Okklusionsaufbau



Telio Activator / Telio Lab Cold Liquid
Activation / Aktivierung

X

X

X

X

X

SR Compositiv
Bonding / Haftvermittler

X

X

X

X

X

Telio Stains (blue, light-yellow, white, orange)
Individualization / Individualisierung

X

Telio Lab LC Incisal
Cut-Back

X

X

Telio Lab LC Base
Relining / Unterfütterung

X

X

Telio Add-On Flow
Modifications / Modifikationen

X

X

SR Gel
Polymerization / Polymerisation

X

X

X

X

X

Telio CS Link
Cementation / Befestigung

X

X

X

X

X